

---

## Immunohistochemistry Protocol for Paraffin-embedded Tissues

---

### 1. Solutions and reagents

- 1.1. Xylene
- 1.2. Ethanol, anhydrous denatured, histological grade (100%, 95%, 70%, 50%)
- 1.3. Washing buffer/TBST: 1X TBS/0.1% Tween-20, pH to 7.6.
- 1.4. Distilled water (dH<sub>2</sub>O)
- 1.5. Antigen Retrieval Solution: 0.01M Sodium Citrate Buffer, pH 6.0

#### To prepare Antigen Retrieval stock solutions:

**10X Stock:** Dissolve 29.4 g sodium citrate trisodium salt dehydrate (C<sub>6</sub>H<sub>5</sub>Na<sub>3</sub>O<sub>7</sub>·2H<sub>2</sub>O) in 1 liter of dH<sub>2</sub>O. Add 5mL Tween-20.

**1X Working Solution:** Mix 200mL 10X stock with 1800mL dH<sub>2</sub>O; pH to 6.0

- 1.6. 3% Hydrogen Peroxide
- 1.7. Blocking Buffer: 10% serum in PBS (serum origin depends on the host of the secondary antibody)
- 1.8. Primary Antibody Diluent: 5% serum in PBS (serum origin depends on the host of the Secondary antibody)
- 1.9. Hematoxylin
- 1.10. Permanent Mounting Medium

---

### 2. Protocol

#### 2.1. Deparaffinization/Rehydration

2.1.1. Heat slides in an oven at 65 °C for 1 hour.

2.1.2. De-paraffinize/hydrate using the following series of washes: two Xylene washes (3 min each), followed by two 100% ethanol rinses (3 min each), followed by 95% ethanol, 70% ethanol, 50% ethanol, 30% ethanol, followed by TBST wash for 3 min on a shaker.

#### 2.2. Antigen Retrieval

*This is recommended Heat Induced Epitope Retrieval (HIER) using Decloaking Chamber/Pressure Cooker. Hot water bath or Microwave with temperature sensor can be also used (protocol would vary depending on the method used).*

2.2.1. Add 500 ml of dH<sub>2</sub>O to Decloaker/Pressure Cooker.

2.2.2. Immerse slides into staining dish containing Antigen Retrieval Solution. Place staining

dish into decloaking chamber.

**2.2.3.** Program to run for 30 seconds at 125° C, followed by 10 seconds at 90° C.

**2.2.4.** Let it cool down to room temperature (10 – 20 minutes).

**2.2.5.** Removes slides and rinse in TBST.

**2.2.6.** Proceed to Staining step.

### **2.3. Staining**

**2.3.1.** Wash slides with TBST for 3 min on a shaker.

**2.3.2.** Inactivate endogenous peroxidase by covering tissue with 3% hydrogen peroxide for 5 min.

**2.3.3.** Wash slides three times with TBST (3 min each on a shaker).

**2.3.4.** Block slides with the blocking solution for 1 hour.

**2.3.5.** Dilute primary antibody in primary antibody diluent per recommendation on data sheet.

**2.3.6.** Apply primary antibody to each section and incubate overnight in the humidified chamber (4 °C).

**2.3.7.** Wash slides three times with TBST (3 min each on a shaker).

**2.3.8.** Apply to each section secondary HRP-conjugated anti-rabbit antibody diluted in the blocking solution per manufacturer's recommendation; incubate for 30 min at room temperature.

**2.3.9.** Wash slides three times with TBST (5 min each on a shaker).

**2.3.10.** Add freshly prepared DAB substrate to the sections and incubate until stain develops (generally 1 min).

**2.3.11.** Rinse sections with water.

**2.3.12.** Counterstain with Hematoxylin (generally 10 seconds).

**2.3.13.** Rinse sections with water.

**2.3.14.** Dehydrate samples using two washes with 100% Ethanol (3 min each), followed by two rinses with Xylene (3 min each).

**2.3.15.** Mount coverslips on slides using permanent mounting medium.